



## **APPROVED CLEANING AND DISINFECTION PRODUCTS**

The following products have been tested and approved by Anacom MedTek for use on pillow speakers and keyboards.

- Mild detergent and warm water
- Denatured Isopropyl alcohol
- PDI super sani-cloth germicidal disposable wipes
- Clorox disinfecting wipes
- Lysol disinfecting wipes
- Spray 9 disinfectant
- Ethylene Oxide Gas @ less than 130 degree Fahrenheit
- Virasept
- Steriplex
- Clorox Pro Quaternary
- Caviwipes
- 3M disinfectant 40L